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Experiment

Capacitors

The charge q on a capacitor's plate is proportional to the potential difference V across the capacitor. We express this relationship with

$$V = \frac{q}{C}$$

where C is a proportionality constant known as the *capacitance*. The unit used for capacitance is the farad, F (1 farad = 1 coulomb/volt).

If a capacitor of capacitance C (in farads), initially charged to a potential V_0 (in volts) is connected across a resistor R (in ohms), a time-dependent current will flow according to Ohm's law. This situation is shown by the RC (resistor-capacitor) circuit below when the switch is turning to 'b' point.


Figure 1

As the charge flows, the charge q on the capacitor is depleted, reducing the potential across the capacitor, which in turn reduces the current. This process creates an exponentially decreasing current, modeled by

$$V(t) = V_0 e^{-t/RC}$$

The rate of the decrease is determined by the product RC , known as the *time constant* of the circuit. A large time constant means that the capacitor will discharge slowly.

In contrast, when the capacitor is charged, the potential across it approaches the final value exponentially, modeled by

$$V(t) = V_0 (1 - e^{-t/RC})$$

The same time constant, RC , describes the rate of charging as well as discharging.

OBJECTIVES

- Experimentally determine the time constant of a resistor-capacitor circuit.

- Compare the time constant to the value predicted from the component values of the resistance and capacitance.
- Measure the potential across a capacitor as a function of time as it discharges.
- Fit an exponential function to the data. One of the fit parameters corresponds to an experimental time constant.

MATERIALS

Mobile device (ipad Mini), Graphical Analysis, Go Direct Voltage sensor
 Connecting wires with clips, Circuit Board with 2200 μF polarized capacitor, 5.6k Ω – 25.6 k Ω resistors, DC power supply, Digital multimeter

PROCEDURE

1. Set up the equipment.
 - a. Connect the circuit as shown in Figure 1.
 - b. Record the values of the resistor and capacitor in your data table, as well as any tolerance values marked on them.
 - c. Launch Graphical Analysis. Connect the Go Direct Voltage Probe to your mobile device.

Connect the clip leads on the voltage probe across the capacitor. **Note:** Connect the red lead to the side of the + terminal, and the black lead to the side of the – terminal respectively.






- d. Adjust the DC power supply about 3.0 V. (Use the indicator of the instrument)
2. Set up Graphical Analysis for triggering and data collection. Click or tap Mode to open Data Collection Settings.
 - a. Change Rate to 10 samples/s and End Collection to 60 s.
 3. Collect data.
 - a. Click or tap Collect to start data collection. (Charge and discharge)
 4. Fit the exponential function, $y = a \cdot \exp(-cx) + b$, to your data.
 - a. Click or tap Graph Tools, , and choose Apply Curve Fit.
 - b. Select Natural Exponent as the curve fit and click or tap Apply.
 - c. Record the value of the fit parameters in your data table. Notice that the c used in the curve fit is not the same as the C representing capacitance. Compare the fit equation to the mathematical model for a capacitor discharge proposed in the introduction.

$$V(t) = V_0 e^{-\frac{t}{RC}}$$

How is fit constant c related to the time constant of the circuit, which was defined in the introduction?

5. Adjust the power supply about 6.0V to 12V and repeat Steps 2–4.

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