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
stopper the separatory funnel and carefully invert it. Without shaking the separatory funnel, immediately open the stopcock to release the pressure. Close the stopcock, shake the funnel several times, and again open the stopcock to release the pressure. Repeat two or three times, venting frequently. Allow the mixture to stand in the separatory funnel until the two layers have completely separated. Carefully drain the lower layer, label as aqueous and save until the final product is isolated.


Extraction/Washing. The operations in this paragraph should be done as rapidly as possible to prevent a reverse reaction and reformation of starting alcohol. Save all washes and extracts until you are sure you have your product. With the organic layer still in the separatory funnel, perform the following manipulations. Wash (twice) and shake the organic layer with 10 mL of water. Drain the lower aqueous layer. Add a 10-mL portion of 10% aqueous sodium bicarbonate to the separatory funnel. Gently swirl the funnel (unscrewed) until the contents are thoroughly mixed observing the formation of carbon dioxide gas. Stopper the funnel, carefully invert, and release the pressure through the stopcock. Gently shake the separatory funnel, with frequent venting for about 1 minute. When minimal pressure buildup is observed, shake the funnel vigorously, venting two or three times. Allow the layers to separate and drain the lower aqueous layer. Wash the organic layer once more with a 10-mL portion of aqueous NaCl (brine) and drain the aqueous layer.

Drying Agent. Pour the organic layer out through the top of the separatory funnel into a small Erlenmeyer flask. Dry the crude product by adding enough anhydrous calcium chloride to lightly cover the bottom of the flask. Gently swirl and allow the drying agent to settle for at least five minutes. Check for signs that the product is dry - the solution is clear, not cloudy, there are no visible droplets of water present, there is not a separate layer of liquid or "puddle", and there is drying agent present that has not "clumped". In other words not needed. If necessary, add more CaCl₂, but be careful not to add too much since product will become absorbed onto the drying agent and reduce the yield.

Distillation. Transfer the clear liquid to a dry 50-mL, round-bottomed flask by decanting or using a Pasteur pipet. Add a boiling stone and distill the crude product using an apparatus as illustrated below, on the technique video, and in Figure 12.7, p. 181 of Mohrig. Use a metal clamp on the distilling flask the technique video, and use the blue Ketch clamps on the ends of the condenser and on the vacuum adapter (ALWAYS) and use the blue Ketch clamps on the ends of the condenser and on the vacuum adapter. Anchor the water out hose in the sink by placing a plastic funnel on the end. Raise the heating mantle off of the benchtop with either a wooden block or ringstand. This precaution allows the heat source to be quickly removed without moving the whole apparatus should something go wrong. Be sure to plug the thermowell heating mantle into the Variac outlet and NOT directly into a regular outlet. Also, never distill to dryness. To prevent loss of product use the vacuum adapter and a pre-weighed 25-mL, round-bottomed flask cooled in an ice bath. Accurately record the boiling point range. The boiling point for this particular compound is usually about 5 °C lower than the literature boiling point. Determine a mass of the product, record its physical characteristics, and calculate a % yield.

Continue to next page.


2311 Simple Distillation Apparatus

Waste Disposal

Dispose of excess 2-butanol in the ORGANIC SOLVENT HAZARDOUS WASTE. Dispose of the pot residue from the distillation into the waste labeled HALOGENATED ORGANIC SOLVENT HAZARDOUS WASTE. Unless otherwise indicated the product will be collected for use in a later reaction. Combine all aqueous washes from the extraction together, test the pH and neutralize with solid baking soda (box) before flushing down the drain with excess water.

Reaction of 2-butanol Alcohol with HCl Jan 26 Tue
Objective: The objective of the experiment is to perform the reaction of 2-butanol with hydrochloric acid and isolate a major liquid product which is predicted to be 2-butyl chloride. The major product is predicted to be a product of substitution and elimination reactions.

Property/Value Table

Property/Value	2-butanol	2-butyl chloride	2-butene	1-butanol
MW	94.12	118.17	98.14	74.12
Color	colorless	colorless	colorless	colorless
solubility	soluble	soluble	soluble	soluble
boiling pt	117.7	117.7	36.9	78.4
density	0.81	0.88	0.68	0.81
refractive index	1.39	1.41	1.39	1.39
miscibility	miscible	miscible	immiscible	miscible
flammable	flammable	flammable	flammable	flammable

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